

For Research Use Only

Mouse Fas/CD95 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84110-2-PBS



Basic Information

Catalog Number: 84110-2-PBS	GenBank Accession Number: NM_007987.2	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 14102	CloneNo.: 241229D3
Source: Rabbit	UNIPROT ID: P25446	
Isotype: IgG	Full Name: Fas (TNF receptor superfamily member 6)	
	Calculated MW: 37kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01025-2, MOUSE FAS/CD95 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84110-1-PBS. Detection antibody: 84110-2-PBS. Standard: Eg1375. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01025-3, MOUSE FAS/CD95 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84110-3-PBS. Detection antibody: 84110-2-PBS. Standard: Eg1375. Range: 1.25-80 ng/mL