

For Research Use Only

IL23R Recombinant antibody, PBS Only (Detector)



Catalog Number: 84112-3-PBS

Basic Information

Catalog Number: 84112-3-PBS	GenBank Accession Number: NM_144701.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 149233	CloneNo.: 240905F2
Source: Rabbit	UNIPROT ID: Q5VWK5-1	
Isotype: IgG	Full Name: interleukin 23 receptor	
	Calculated MW: 82 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01048-3, IL23R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84112-4-PBS. Detection antibody: 84112-3-PBS. Standard: Eg0853. Range: 0.625-40 ng/mL