

For Research Use Only

# Rat gp130/IL6ST Recombinant antibody, PBS Only (Capture)



Catalog Number: 84113-4-PBS

## Basic Information

Catalog Number: 84113-4-PBS	GenBank Accession Number: NM_001008725.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 25205	CloneNo.: 241340G7
Source: Rabbit	Full Name: interleukin 6 signal transducer	
Isotype: IgG	Calculated MW: 102 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
rat

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01049-2, RAT gp130 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84113-4-PBS. Detection antibody: 84113-1-PBS. Standard: Eg0857. Range: 0.625-80 ng/mL.