

For Research Use Only

PPID Recombinant antibody, PBS Only (Capture)

Catalog Number: 84115-1-PBS



Basic Information

Catalog Number: 84115-1-PBS	GenBank Accession Number: BC030707	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5481	CloneNo.: 241361B2
Source: Rabbit	UNIPROT ID: Q08752	
Isotype: IgG	Full Name: peptidylprolyl isomerase D	
Immunogen Catalog Number: AG3409	Calculated MW: 370 aa, 41 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01018-1, PPID Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84115-1-PBS. Detection antibody: 84115-5-PBS. Standard: Ag3409. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01018-2, PPID Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84115-1-PBS. Detection antibody: 84115-2-PBS. Standard: Ag3409. Range: 0.625-80 ng/mL