

For Research Use Only

CHRN2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84116-4-PBS

Basic Information

Catalog Number: 84116-4-PBS	GenBank Accession Number: BC075041	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1141	CloneNo.: 241374G9
Source: Rabbit	UNIPROT ID: P17787	
Isotype: IgG	Full Name: cholinergic receptor, nicotinic, beta 2 (neuronal)	
Immunogen Catalog Number: AG12208	Calculated MW: 502 aa, 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01019-2, CHRNB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84116-1-PBS. Detection antibody: 84116-4-PBS. Standard: Ag12208. Range: 0.625-80 ng/mL.