

For Research Use Only

# ADAM20 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84122-2-PBS

## Basic Information

Catalog Number: 84122-2-PBS	GenBank Accession Number: BC025378	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8748	CloneNo.: 241149H7
Source: Rabbit	UNIPROT ID: O43506	
Isotype: IgG	Full Name: ADAM metallopeptidase domain 20	
Immunogen Catalog Number: AG25814	Calculated MW: 776 aa, 87 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01020-1, ADAM20 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84122-2-PBS. Detection antibody: 84122-1-PBS. Standard: Ag25814. Range: 0.781-100 ng/mL