

For Research Use Only

HTR2A Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 84126-1-PBS

Basic Information

Catalog Number: 84126-1-PBS	GenBank Accession Number: BC074848	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3356	CloneNo.: 241440C12
Source: Rabbit	UNIPROT ID: P28223	
Isotype: IgG	Full Name: 5-hydroxytryptamine (serotonin) receptor 2A	
Immunogen Catalog Number: AG24056	Calculated MW: 471 aa, 53 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01029-2, 5HT2A Receptor Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84126-1-PBS. Detection antibody: 84126-3-PBS. Standard: Ag24056. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP01029-3, 5HT2A Receptor Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84126-2-PBS. Detection antibody: 84126-1-PBS. Standard: Ag24056. Range: 0.625-40 ng/mL.