

For Research Use Only

DPP10 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84127-1-PBS

Basic Information

Catalog Number: 84127-1-PBS	GenBank Accession Number: BC030832	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57628	CloneNo.: 241446A1
Source: Rabbit	UNIPROT ID: Q8N608	
Isotype: IgG	Full Name: dipeptidyl-peptidase 10	
Immunogen Catalog Number: AG3231	Calculated MW: 796 aa, 91 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01035-2, DPP10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84127-1-PBS. Detection antibody: 84127-4-PBS. Standard: Ag3231. Range: 0.625-80 ng/mL