

For Research Use Only

# DPP10 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84127-4-PBS

## Basic Information

Catalog Number: 84127-4-PBS	GenBank Accession Number: BC030832	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57628	CloneNo.: 241446H3
Source: Rabbit	UNIPROT ID: Q8N608	
Isotype: IgG	Full Name: dipeptidyl-peptidase 10	
Immunogen Catalog Number: AG3231	Calculated MW: 796 aa, 91 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01035-2, DPP10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84127-1-PBS. Detection antibody: 84127-4-PBS. Standard: Ag3231. Range: 0.625-80 ng/mL