

For Research Use Only

RAB9A Recombinant antibody, PBS Only (Detector)



Catalog Number: 84129-1-PBS

Basic Information

Catalog Number: 84129-1-PBS	GenBank Accession Number: BC017265	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9367	CloneNo.: 241372A6
Source: Rabbit	UNIPROT ID: P51151	
Isotype: IgG	Full Name: RAB9A, member RAS oncogene family	
Immunogen Catalog Number: AG1969	Calculated MW: 201 aa, 23 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01043-3, RAB9A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84129-4-PBS. Detection antibody: 84129-1-PBS. Standard: Ag1969. Range: 0.313-40 ng/mL