

For Research Use Only

TMIGD1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84130-4-PBS



Basic Information

Catalog Number: 84130-4-PBS	GenBank Accession Number: BC137201	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 388364	CloneNo.: 241399G9
Source: Rabbit	UNIPROT ID: Q6UXZ0	
Isotype: IgG	Full Name: transmembrane and immunoglobulin domain containing 1	
Immunogen Catalog Number: AG25691	Calculated MW: 262 aa, 29 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP01037-1, TMIGD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84130-4-PBS. Detection antibody: 84130-3-PBS. Standard: Ag25691. Range: 1.563-200 ng/mL



Cytometric bead array standard curve of MP01037-2, TMIGD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84130-1-PBS. Detection antibody: 84130-4-PBS. Standard: Ag25691. Range: 1.563-200 ng/mL



Biolayer interferometry (BLI) kinetic assay of 84130-4-PBS against Human TMIGD1 was performed. The affinity constant is 4.19 nM.