

For Research Use Only

# TTC18 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84133-4-PBS

## Basic Information

Catalog Number: 84133-4-PBS	GenBank Accession Number: BC032856	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 118491	CloneNo.: 241177E11
Source: Rabbit	UNIPROT ID: Q5TON1	
Isotype: IgG	Full Name: tetra(ri)copeptide repeat domain 18	
Immunogen Catalog Number: AG35209	Calculated MW: 126 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01054-1, TTC18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84133-1-PBS. Detection antibody: 84133-4-PBS. Standard: Ag35209. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01054-3, TTC18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84133-3-PBS. Detection antibody: 84133-4-PBS. Standard: Ag35209. Range: 0.781-100 ng/mL.