

For Research Use Only

TTC18 Recombinant antibody

Catalog Number: 84133-5-RR



Basic Information

Catalog Number: 84133-5-RR	GenBank Accession Number: BC032856	Purification Method: Protein A purification
Size: 1000 ug/ml	GeneID (NCBI): 118491	CloneNo.: 241177C9
Source: Rabbit	UNIPROT ID: Q5TON1	Recommended Dilutions: WB 1:1000-1:4000
Isotype: IgG	Full Name: tetra(ri)copeptide repeat domain 18	
Immunogen Catalog Number: AG35209	Calculated MW: 126 kDa	
	Observed MW: 126 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HaCaT cells,
Species Specificity: human	

Background Information

Cilia- and flagella-associated protein 70 (CFAP70) harbours eight tetra(ri)copeptide repeats (TPRs) and bound tightly to the ciliary axoneme. CFAP70 is necessary to assemble spermatid flagella and could be used as a diagnostic target for male infertility with OAT in the clinic. CFAP70 is a novel regulatory component of the ODA in motile cilia and flagella, and that the N-terminus is important for its ciliary localization (PMID: 30158508).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84133-5-RR (TTC18 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 84133-5-RR against Human TTC18 were performed. The affinity constant is 5.06 nM.