

For Research Use Only

# PCDH8 Recombinant antibody, PBS Only (Capture)



Catalog Number: 84134-4-PBS

## Basic Information

|                                     |                                       |  |
|-------------------------------------|---------------------------------------|--|
| Catalog Number:<br>84134-4-PBS      | GenBank Accession Number:<br>BC036025 | Purification Method:<br>Protein A purification |
| Size:<br>1 mg/ml                    | GeneID (NCBI):<br>5100                | CloneNo.:<br>241362F11                         |
| Source:<br>Rabbit                   | UNIPROT ID:<br>O95206                 |  |
| Isotype:<br>IgG                     | Full Name:<br>protocadherin 8         |  |
| Immunogen Catalog Number:<br>AG5090 | Calculated MW:<br>1070 aa, 113 kDa    |  |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01033-1, PCDH8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84134-4-PBS. Detection antibody: 84134-3-PBS. Standard: Ag5090. Range: 0.625-80 ng/mL