

For Research Use Only

NELF-B Recombinant antibody, PBS Only (Detector)



Catalog Number: 84135-1-PBS

Basic Information

Catalog Number: 84135-1-PBS	GenBank Accession Number: BC011892	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 25920	CloneNo.: 24143486
Source: Rabbit	UNIPROT ID: Q8WX92	
Isotype: IgG	Full Name: cofactor of BRCA 1	
Immunogen Catalog Number: AG9599	Calculated MW: 580aa,66 kDa; 475aa,54 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01055-1, NELF-B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84135-4-PBS. Detection antibody: 84135-1-PBS. Standard: Ag9599. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01055-2, NELF-B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84135-2-PBS. Detection antibody: 84135-1-PBS. Standard: Ag9599. Range: 0.781-100 ng/mL.