

For Research Use Only

Mouse IL-1 beta/IL-1F2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84136-4-PBS

Basic Information

Catalog Number:

84136-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_008361.4

GeneID (NCBI):

16176

UNIPROT ID:

P10749

Full Name:

interleukin 1 beta

Calculated MW:

31kDa

Purification Method:

Protein A purification

CloneNo.:

241364E12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01046-2, MOUSE IL-1 beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84136-2-PBS. Detection antibody: 84136-4-PBS. Standard: Eg1577. Range: 0.313-10 ng/mL.