

For Research Use Only

Mouse CD80/B7-1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84137-2-PBS



Basic Information

Catalog Number: 84137-2-PBS	GenBank Accession Number: NM_009855.2	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 12519	CloneNo.: 241371C5
Source: Rabbit	UNIPROT ID: Q00609-1	
Isotype: IgG	Full Name: CD80 antigen	
	Calculated MW: 35 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01051-1, MOUSE CD80 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84137-2-PBS. Detection antibody: 84137-4-PBS. Standard: Eg0840. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01051-3, MOUSE CD80 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84137-2-PBS. Detection antibody: 84137-3-PBS. Standard: Eg0840. Range: 0.313-40 ng/mL