

For Research Use Only

CD155/PVR Recombinant antibody, PBS Only (Detector)

Catalog Number: 84138-4-PBS



Basic Information

Catalog Number: 84138-4-PBS	GenBank Accession Number: BC015542	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5817	CloneNo.: 241405E7
Source: Rabbit	UNIPROT ID: P15151	
Isotype: IgG	Full Name: poliovirus receptor	
	Calculated MW: 45 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01053-2, CD155/PVR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84138-1-PBS. Detection antibody: 84138-4-PBS. Standard: Eg0971. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP01053-3, CD155/PVR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84138-3-PBS. Detection antibody: 84138-4-PBS. Standard: Eg0971. Range: 0.313-40 ng/mL.