

For Research Use Only

CD3EAP Recombinant antibody, PBS Only (Detector)



Catalog Number: 84139-3-PBS

Basic Information

Catalog Number: 84139-3-PBS	GenBank Accession Number: BC038992	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10849	CloneNo.: 241379E5
Source: Rabbit	UNIPROT ID: O15446	
Isotype: IgG	Full Name: CD3e molecule, epsilon associated protein	
Immunogen Catalog Number: AG35074	Calculated MW: 510 aa, 55 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01056-1, CD3EAP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84139-2-PBS. Detection antibody: 84139-3-PBS. Standard: Ag35074. Range: 0.625-80 ng/mL.