

For Research Use Only

# AMH Recombinant antibody, PBS Only (Detector)

Catalog Number: 84153-1-PBS



## Basic Information

Catalog Number: 84153-1-PBS	GenBank Accession Number: BC049194	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 268	CloneNo.: 241148F7
Source: Rabbit	UNIPROT ID: P03971	
Isotype: IgG	Full Name: anti-Mullerian hormone	
Immunogen Catalog Number: AG5670	Calculated MW: 59 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01075-1, AMH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84153-2-PBS. Detection antibody: 84153-1-PBS. Standard: Ag5670. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01075-3, AMH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84153-4-PBS. Detection antibody: 84153-1-PBS. Standard: Ag5670. Range: 0.781-100 ng/mL.