

For Research Use Only

Mouse CD27 Recombinant antibody, PBS Only (Capture)

Catalog Number:84157-1-PBS



Basic Information

Catalog Number:

84157-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001033126.2

GeneID (NCBI):

21940

UNIPROT ID:

P41272

Full Name:

CD27 antigen

Calculated MW:

28kDa

Purification Method:

Protein A purification

CloneNo.:

241378A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01069-3, MOUSE CD27/TNFRSF7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84157-1-PBS. Detection antibody: 84157-4-PBS. Standard: Eg1311. Range: 0.781-50 ng/mL