

For Research Use Only

ARL4D Recombinant antibody, PBS Only (Capture)

Catalog Number: 84163-2-PBS



Basic Information

| | | |
|--|--|---|
| Catalog Number: 84163-2-PBS | GenBank Accession Number: BC000043 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 379 | CloneNo.: 241503E12 |
| Source: Rabbit | UNIPROT ID: P49703 | |
| Isotype: IgG | Full Name: ADP-ribosylation factor-like 4D | |
| Immunogen Catalog Number: AG7369 | Calculated MW: 22 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01065-1, ARL4D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84163-2-PBS. Detection antibody: 84163-1-PBS. Standard: Ag7369. Range: 1.563-200 ng/mL.