

For Research Use Only

ABCC10 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84164-1-PBS



Basic Information

Catalog Number: 84164-1-PBS	GenBank Accession Number: NM_001198934.1	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 89845	CloneNo.: 241429C7
Source: Rabbit	UNIPROT ID: Q5T3U5	
Isotype: IgG	Full Name: ATP-binding cassette, sub-family C (CFTR/MRP), member 10	
Immunogen Catalog Number: AG33719	Calculated MW: 161kd	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01067-1, ABCC10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84164-1-PBS. Detection antibody: 84164-3-PBS. Standard: Ag33719. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 84164-1-PBS against Human ABCC10 was performed. The affinity constant is below 1 pM.