

For Research Use Only

# IL-7 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84165-4-PBS



## Basic Information

Catalog Number:

84165-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_000880.4

GeneID (NCBI):

3574

UNIPROT ID:

P13232-1

Full Name:

interleukin 7

Calculated MW:

20 kDa

Purification Method:

Protein A purification

CloneNo.:

241000G11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01089-1, IL-7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84165-4-PBS. Detection antibody: 84165-2-PBS. Standard: Eg0830. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP01089-3, IL-7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84165-4-PBS. Detection antibody: 84165-3-PBS. Standard: Eg0830. Range: 0.313-20 ng/mL