

For Research Use Only

PCSK9 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84172-1-PBS



Basic Information

Catalog Number: 84172-1-PBS	GenBank Accession Number: NM_174936.4	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 255738	CloneNo.: 241395B8
Source: Rabbit	UNIPROT ID: Q8NBP7	
Isotype: IgG	Full Name: proprotein convertase subtilisin/kexin type 9	
	Calculated MW: 74 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01071-1, PCSK9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84172-1-PBS. Detection antibody: 84172-5-PBS. Standard: Eg0362. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP01071-2, PCSK9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84172-1-PBS. Detection antibody: 84172-2-PBS. Standard: Eg0362. Range: 0.781-50 ng/mL