

For Research Use Only

# DPP3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84173-1-PBS



## Basic Information

Catalog Number: 84173-1-PBS	GenBank Accession Number: BC007221	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10072	CloneNo.: 241360B7
Source: Rabbit	UNIPROT ID: Q9NY33	
Isotype: IgG	Full Name: dipeptidyl-peptidase 3	
Immunogen Catalog Number: AG1022	Calculated MW: 83 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01074-1, DPP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84173-1-PBS. Detection antibody: 84173-2-PBS. Standard: Ag1022. Range: 1.563-200 ng/mL.