

For Research Use Only

# Kininogen-1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84175-1-PBS



## Basic Information

Catalog Number: 84175-1-PBS	GenBank Accession Number: NM_001102416.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3827	CloneNo.: 241397A7
Source: Rabbit	UNIPROT ID: P01042-1	
Isotype: IgG	Full Name: kininogen 1	
	Calculated MW: 72kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01081-3, Kininogen-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84175-4-PBS. Detection antibody: 84175-1-PBS. Standard: Eg1162. Range: 0.781-100 ng/mL.