

For Research Use Only

# ERBB2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84182-1-PBS



## Basic Information

Catalog Number:

84182-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC156755

GeneID (NCBI):

2064

UNIPROT ID:

P04626

Full Name:

v-erb-b2 erythroblastic leukemia  
viral oncogene homolog 2,  
neuro/glioblastoma derived  
oncogene homolog (avian)

Calculated MW:

1255 aa, 138 kDa

Purification Method:

Protein A purification

CloneNo.:

241409D2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01080-2, ERBB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84182-2-PBS. Detection antibody: 84182-1-PBS. Standard: RPO2350. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01080-3, ERBB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84182-3-PBS. Detection antibody: 84182-1-PBS. Standard: RPO2350. Range: 0.625-80 ng/mL