

For Research Use Only

4-1BBL/TNFSF9 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84185-1-PBS



Basic Information

Catalog Number:

84185-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC104805

GeneID (NCBI):

8744

UNIPROT ID:

P41273

Full Name:

tumor necrosis factor (ligand)
superfamily, member 9

Calculated MW:

254 aa, 27 kDa

Purification Method:

Protein A purification

CloneNo.:

241386B8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01091-2, 4-1BBL/TNFSF9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84185-4-PBS. Detection antibody: 84185-1-PBS. Standard: Eg0010. Range: 0.313-20 ng/mL.