

For Research Use Only

# 4-1BBL/TNFSF9 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84185-3-PBS



## Basic Information

Catalog Number:

84185-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC104805

GeneID (NCBI):

8744

UNIPROT ID:

P41273

Full Name:

tumor necrosis factor (ligand)  
superfamily, member 9

Calculated MW:

254 aa, 27 kDa

Purification Method:

Protein A purification

CloneNo.:

241386D5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01091-1, 4-1BBL/TNFSF9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84185-2-PBS. Detection antibody: 84185-3-PBS. Standard: Eg0010. Range: 0.156-20 ng/mL.