

For Research Use Only

# t-Plasminogen activator/tPA Recombinant antibody, PBS Only (Capture)

Catalog Number: 84191-1-PBS



## Basic Information

Catalog Number:	84191-1-PBS	GenBank Accession Number:	BC095403	Purification Method:	Protein A purification
Size:	1 mg/ml	GeneID (NCBI):	5327	CloneNo.:	241526C5
Source:	Rabbit	UNIPROT ID:	P00750		
Isotype:	IgG	Full Name:	plasminogen activator, tissue		

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01082-1, PLAT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84191-1-PBS. Detection antibody: 84191-4-PBS. Standard: Eg0178. Range: 0.781-100 ng/mL.