

For Research Use Only

# LIMPII Recombinant antibody

Catalog Number: 84193-4-RR



## Basic Information

Catalog Number:

84193-4-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25908

GenBank Accession Number:

BC021892

GeneID (NCBI):

950

UNIPROT ID:

Q14108

Full Name:

scavenger receptor class B, member 2

Calculated MW:

478 aa, 54 kDa

Observed MW:

75 kDa

Purification Method:

Protein A purification

CloneNo.:

241468F9

Recommended Dilutions:

WB 1:2000-1:10000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, MCF-7 cells, Y79 cells, SH-SY5Y cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84193-4-RR (LIMPII antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 84193-4-RR against Human LIMPII were performed. The affinity constant is below 1 pM.