

For Research Use Only

ATAD1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84194-2-PBS



Basic Information

Catalog Number: 84194-2-PBS	GenBank Accession Number: BC010868	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 84896	CloneNo.: 241449C7
Source: Rabbit	UNIPROT ID: Q8NBU5	
Isotype: IgG	Full Name: ATPase family, AAA domain containing 1	
Immunogen Catalog Number: AG10343	Calculated MW: 361 aa, 41 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP01090-1, ATAD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84194-4-PBS. Detection antibody: 84194-2-PBS. Standard: Ag10343. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01090-2, ATAD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84194-3-PBS. Detection antibody: 84194-2-PBS. Standard: Ag10343. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01090-3, ATAD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84194-1-PBS. Detection antibody: 84194-2-PBS. Standard: Ag10343. Range: 0.625-80 ng/mL