

For Research Use Only

METTL6 Recombinant antibody, PBS Only (Capture)

Catalog Number:84209-3-PBS



Basic Information

Catalog Number: 84209-3-PBS	GenBank Accession Number: BC022400	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 131965	CloneNo.: 241482G12
Source: Rabbit	UNIPROT ID: Q8TCB7	
Isotype: IgG	Full Name: methyltransferase like 6	
Immunogen Catalog Number: AG9714	Calculated MW: 255 aa, 30 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01116-1, METTL6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84209-3-PBS. Detection antibody: 84209-2-PBS. Standard: Ag9714. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01116-3, METTL6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84209-3-PBS. Detection antibody: 84209-1-PBS. Standard: Ag9714. Range: 1.25-40 ng/mL