

For Research Use Only

ANKRD6 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84210-2-PBS



Basic Information

Catalog Number: 84210-2-PBS	GenBank Accession Number: BC042173	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 22881	CloneNo.: 241494H2
Source: Rabbit	UNIPROT ID: Q9Y2G4	
Isotype: IgG	Full Name: ankyrin repeat domain 6	
Immunogen Catalog Number: AG19328	Calculated MW: 727 aa, 80 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01124-1, ANKRD6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84210-2-PBS. Detection antibody: 84210-1-PBS. Standard: Ag19328. Range: 0.781-100 ng/mL