

For Research Use Only

ANKRD6 Recombinant antibody, PBS Only

Catalog Number: 84210-4-PBS



Basic Information

Catalog Number: 84210-4-PBS	GenBank Accession Number: BC042173	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 22881	CloneNo.: 241494D11
Source: Rabbit	UNIPROT ID: Q9Y2G4	
Isotype: IgG	Full Name: ankyrin repeat domain 6	
Immunogen Catalog Number: AG19328	Calculated MW: 727 aa, 80 kDa	

Applications

Tested Applications:
IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed A549 cells using ANKRD6 antibody (84210-4-RR, Clone: 241494D11) at dilution of 1:250 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 84210-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84210-4-RR against Human ANKRD6 were performed. The affinity constant is 0.722 nM.