

For Research Use Only

ANKRD6 Recombinant antibody

Catalog Number: 84210-4-RR



Basic Information

Catalog Number: 84210-4-RR	GenBank Accession Number: BC042173	Purification Method: Protein A purification
Size: 900 ug/ml	GeneID (NCBI): 22881	CloneNo.: 241494D11
Source: Rabbit	UNIPROT ID: Q9Y2G4	Recommended Dilutions: IF/ICC 1:125-1:500
Isotype: IgG	Full Name: ankyrin repeat domain 6	
Immunogen Catalog Number: AG19328	Calculated MW: 727 aa, 80 kDa	

Applications

Tested Applications: IF/ICC, ELISA	Positive Controls: IF/ICC : A549 cells,
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed A549 cells using ANKRD6 antibody (84210-4-RR, Clone: 241494D11) at dilution of 1:250 and CoraLite@488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).



Biolayer interferometry (BLI) kinetic assays of 84210-4-RR against Human ANKRD6 were performed. The affinity constant is 0.722 nM.