

For Research Use Only

# SNX18 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84215-3-PBS



## Basic Information

|                                      |                                       |  |
|--------------------------------------|---------------------------------------|--|
| Catalog Number:<br>84215-3-PBS       | GenBank Accession Number:<br>BC060791 | Purification Method:<br>Protein A purification |
| Size:<br>1 mg/ml                     | GeneID (NCBI):<br>112574              | CloneNo.:<br>241537C6                          |
| Source:<br>Rabbit                    | UNIPROT ID:<br>Q96RF0                 |  |
| Isotype:<br>IgG                      | Full Name:<br>sorting nexin 18        |  |
| Immunogen Catalog Number:<br>AG16645 | Calculated MW:<br>628 aa, 69 kDa      |  |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01150-1, SNX18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84215-1-PBS. Detection antibody: 84215-3-PBS. Standard: Ag16645. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01150-2, SNX18 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84215-2-PBS. Detection antibody: 84215-3-PBS. Standard: Ag16645. Range: 1.25-40 ng/mL