

For Research Use Only

KLHDC3 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84217-2-PBS



Basic Information

Catalog Number: 84217-2-PBS	GenBank Accession Number: BC041793	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 116138	CloneNo.: 241475B12
Source: Rabbit	UNIPROT ID: Q9BQ90	
Isotype: IgG	Full Name: kelch domain containing 3	
Immunogen Catalog Number: AG4879	Calculated MW: 382 aa, 43 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01151-1, KLHDC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84217-2-PBS. Detection antibody: 84217-3-PBS. Standard: Ag4879. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01151-2, KLHDC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84217-4-PBS. Detection antibody: 84217-2-PBS. Standard: Ag4879. Range: 0.781-100 ng/mL