

For Research Use Only

FSP1/S100A4 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number:84221-2-PBS



Basic Information

Catalog Number: 84221-2-PBS	GenBank Accession Number: BC016300	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6275	CloneNo.: 241448B1
Source: Rabbit	ENSEMBL Gene ID: ENSG00000196154	
Isotype: IgG	UNIPROT ID: P26447	
Immunogen Catalog Number: AG9019	Full Name: S100 calcium binding protein A4	
	Calculated MW: 101 aa, 12 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP01155-1, S100A4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84221-3-PBS. Detection antibody: 84221-2-PBS. Standard: Ag9019. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01155-2, S100A4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84221-2-PBS. Detection antibody: 84221-1-PBS. Standard: Ag9019. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01155-3, S100A4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84221-2-PBS. Detection antibody: 84221-4-PBS. Standard: Ag9019. Range: 0.625-80 ng/mL