

For Research Use Only

C19orf12 Recombinant antibody, PBS Only (Detector)

Catalog Number:84223-1-PBS



Basic Information

Catalog Number: 84223-1-PBS	GenBank Accession Number: BC063518	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 83636	CloneNo.: 241544G9
Source: Rabbit	UNIPROT ID: Q9NSK7	
Isotype: IgG	Full Name: chromosome 19 open reading frame 12	
Immunogen Catalog Number: AG18308	Calculated MW: 107 aa, 11 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01161-1, C19orf12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84223-2-PBS. Detection antibody: 84223-1-PBS. Standard: Ag18308. Range: 1.563-100 ng/mL.