

For Research Use Only

C19orf12 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84223-2-PBS



Basic Information

Catalog Number:

84223-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18308

GenBank Accession Number:

BC063518

GeneID (NCBI):

83636

UNIPROT ID:

Q9NSK7

Full Name:

chromosome 19 open reading frame

12

Calculated MW:

107 aa, 11 kDa

Purification Method:

Protein A purification

CloneNo.:

241544G11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01161-1, C19orf12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84223-2-PBS. Detection antibody: 84223-1-PBS. Standard: Ag18308. Range: 1.563-100 ng/mL.