

For Research Use Only

# ZC3H13 Recombinant antibody

Catalog Number: 84228-4-RR



## Basic Information

<b>Catalog Number:</b> 84228-4-RR	<b>GenBank Accession Number:</b> NM_001076788	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1000 µg/ml	<b>GeneID (NCBI):</b> 23091	<b>CloneNo.:</b> 241516G3
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q5T200	<b>Recommended Dilutions:</b> WB 1:500-1:2000
<b>Isotype:</b> IgG	<b>Full Name:</b> zinc finger CCCH-type containing 13	
<b>Immunogen Catalog Number:</b> AG28909	<b>Calculated MW:</b> 197 kDa	
	<b>Observed MW:</b> 196-240 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : NIH/3T3 cells, HeLa cells
<b>Species Specificity:</b> human, mouse	

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84228-4-RR (ZC3H13 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.