

For Research Use Only

# CD122/IL-2RB Recombinant antibody, PBS Only (Capture)

Catalog Number:84229-3-PBS



## Basic Information

<b>Catalog Number:</b> 84229-3-PBS	<b>GenBank Accession Number:</b> BC025691	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 3560	<b>CloneNo.:</b> 241408B1
<b>Source:</b> Rabbit	<b>ENSEMBL Gene ID:</b> ENSG00000100385	
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> P14784	
<b>Immunogen Catalog Number:</b> AG3929	<b>Full Name:</b> interleukin 2 receptor, beta	
	<b>Calculated MW:</b> 551 aa, 61 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01123-1, CD122/IL2RB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84229-3-PBS. Detection antibody: 84229-2-PBS. Standard: Ag3929. Range: 0156-20 ng/mL.



Cytometric bead array standard curve of MP01123-2, CD122/IL2RB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84229-3-PBS. Detection antibody: 84229-1-PBS. Standard: Ag3929. Range: 0156-20 ng/mL.