

For Research Use Only

TNIK Recombinant antibody, PBS Only (Detector)

Catalog Number: 84231-3-PBS



Basic Information

Catalog Number: 84231-3-PBS	GenBank Accession Number: BC150256	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 23043	CloneNo.: 241518H3
Source: Rabbit	UNIPROT ID: Q9UKE5	
Isotype: IgG	Full Name: TRAF2 and NCK interacting kinase	
Immunogen Catalog Number: AG34132		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01127-2, TNK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84231-2-PBS. Detection antibody: 84231-3-PBS. Standard: Ag34132. Range: 0.781-100 ng/mL.