

For Research Use Only

LIF Recombinant antibody, PBS Only (Capture)

Catalog Number: 84232-3-PBS



Basic Information

| | | |
|---------------------------------------|---|---|
| Catalog Number: 84232-3-PBS | GenBank Accession Number: BC093733 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 3976 | CloneNo.: 241519A7 |
| Source: Rabbit | UNIPROT ID: P15018 | |
| Isotype: IgG | Full Name: leukemia inhibitory factor (cholinergic differentiation factor) | |
| | Calculated MW: 202 aa, 22 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01130-1, LIF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84232-3-PBS. Detection antibody: 84232-2-PBS. Standard: Eg0111. Range: 0.625-80 ng/mL