

For Research Use Only

# CD79B Recombinant antibody, PBS Only (Capture)



Catalog Number: 84234-3-PBS

## Basic Information

Catalog Number:

84234-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC032651

GeneID (NCBI):

974

UNIPROT ID:

P40259

Full Name:

CD79b molecule, immunoglobulin-associated beta

Calculated MW:

229 aa, 26 kDa

Purification Method:

Protein A purification

CloneNo.:

241274C6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01141-1, CD79B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84234-3-PBS. Detection antibody: 84234-2-PBS. Standard: RP02310. Range: 0.156-20 ng/mL