

For Research Use Only

APOBEC3A Recombinant antibody, PBS Only (Detector)

Catalog Number: 84239-5-PBS



Basic Information

Catalog Number: 84239-5-PBS	GenBank Accession Number: BC126416	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 200315	CloneNo.: 241529H5
Source: Rabbit	UNIPROT ID: P31941	
Isotype: IgG	Full Name: apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3A	
Immunogen Catalog Number: AG18845	Calculated MW: 199 aa, 23 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01122-1, APOBEC3A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84239-1-PBS. Detection antibody: 84239-5-PBS. Standard: Ag18845. Range: 0.625-80 ng/mL.