

For Research Use Only

CD34 Recombinant antibody, PBS Only

Catalog Number: 84242-5-PBS



Basic Information

| | | |
|---------------------------------------|--|---|
| Catalog Number: 84242-5-PBS | GenBank Accession Number: BC039146 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 947 | CloneNo.: 242069F8 |
| Source: Rabbit | UNIPROT ID: P28906 | |
| Isotype: IgG | Full Name: CD34 molecule | |
| | Calculated MW: 41 kDa | |
| | Observed MW: 105 kDa | |

Applications

Tested Applications:
WB, IHC, Indirect ELISA

Species Specificity:
human

Background Information

CD34 is a 105- to 120-kDa glycoposphoprotein expressed on the majority of hematopoietic stem/progenitor cells, bone marrow stromal cells, capillary endothelial cells, embryonic fibroblasts, and some nerve tissue. CD34 is a commonly used marker for identifying human hematopoietic stem/progenitor cells and mediates cell adhesion and lymphocyte homing by binding L-selectin and E-selectin ligands. CD34 is also one of the best negative selection markers for characterizing and/or isolating human MSCs from bone marrow and other sources. Along with other positive selection markers (such as CD29, CD44, CD90, CD105 and CD166), negative selection markers (such as CD34 and CD45) are used for MSC identification.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 84242-5-RR (CD34 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84242-5-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using 84242-5-RR (CD34 antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 84242-5-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using 84242-5-RR (CD34 antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 84242-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84242-5-RR against Human CD34 were performed. The affinity constant is 16.3 nM.