

For Research Use Only

# IFN-gamma Recombinant antibody, PBS Only (Capture)

Catalog Number:84246-1-PBS



## Basic Information

<b>Catalog Number:</b> 84246-1-PBS	<b>GenBank Accession Number:</b> BC070256	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 3458	<b>CloneNo.:</b> 241390A9
<b>Source:</b> Rabbit	<b>ENSEMBL Gene ID:</b> ENSG00000111537	
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> P01579	
	<b>Full Name:</b> IFN gamma	
	<b>Calculated MW:</b> 166 aa, 19 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01131-2, IFNG/IFN-gamma Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84246-1-PBS. Detection antibody: 84246-2-PBS. Standard: Eg0004. Range: 0.625-40 ng/mL